



3-MODE CONTINUOUS

The SMP 3-Mode Continuous Vacuum Regulator is a perfect combination of clearly marked indicators and precise calibration, providing the user with a high quality product for maximum performance and minimum maintenance. In addition, this vacuum regulator provides a broad selection of fittings ensuring compatibility with existing vacuum outlets in any medical facility.



Glow in the dark gauge



Optional high impact safety trap



Color coded flow range

FEATURES AND BENEFITS

Gauge accuracy: ±1.6%

Color coded flow ranges

Glow in the dark gauge for easy reading under low light conditions

Gauge available in ISO and US configurations

Accurate, precise suction control

Ergonomic design

Can be sterilized internally with cold disinfectant

Durable high quality ABS protective body for high impact resistance

Comes pre-equipped with an overflow safety trap made from high performance precision thermoplastic

Available gauges:

(0 - 160 mmHg) / (0 - 200 mmHg) / (0 - 300 mmHg) / (0 - 400 mmHg) / (0 - 760 mmHg)


Shell Color: White or Yellow

Industry leading 7 year warranty



OFF	\checkmark	\checkmark	\checkmark	\checkmark
REGULATE	\checkmark	\checkmark	\checkmark	\checkmark
FULL	\checkmark	\checkmark	\checkmark	\checkmark
RANGE	0 - Full Line Vacuum			
GAUGE DISPLAY RANGE	0 - 200 mmHg	0 - 300 mmHg	0 - 400 mmHg	0 - 760 mmHg
COLOR CODED FOR PRESSURE LEVEL	✓	\checkmark	✓	✓
GAUGE ACCURACY	±1.6%	±1.6%	±1.6%	±1.6%
GAUGE DIRECTION	Clockwise / Counter Clockwise	Clockwise / Counter Clockwise	Counter Clockwise	Clockwise / Counter Clockwise
FLOW RATE	0 - 60L / min without fitting at full pressure	0 - 60L / min without fitting at full pressure	0 - 60L / min without fitting at full pressure	0 - 60L / min without fitting at full pressure
MEDICAL CUICTION				
MEDICAL SUCTION EQUIPMENT STANDARD	ISO 10079-3	ISO 10079-3	ISO 10079-3	ISO 10079-3
	ISO 10079-3 7 Years	ISO 10079-3 7 Years	ISO 10079-3 7 Years	ISO 10079-3 7 Years

ORDERING GUIDE



(PB = PURITAN-BENNETT)



Example:

VYC - 300 - U - DS - T

Color: Yellow shell color Mode: CONTINUOUS Gauge Range: 0 - 300 mmHg Gauge Direction: Clockwise Adapter: DISS adapter

Fitting: Trap

* 0-200 mmHg available in clockwise gauge only

* 0-400 mmHg available in
counter clockwise gauge only

* Intermittent vacuum regulatory

available in white shell only

MODE SETTINGS

REG:

Allows degree of suction to be adjusted by use of the regulating knob



Suction is no longer on or being supplied to the patient



Maximum suction is administered to the patient

SMP