

SMP-HC1000







Full nurse control panel



BED POSITION













DED POSITION

FEATURES AND BENEFITS

Ultra-low height (250 mm)

4 LINAK motors

Power operated back-rest, leg-rest and height control

Power operated trendelenburg, reverse trendelenburg, fowler and vascular position

Collapsible aluminum side rails

Premium HPL head / foot boards and bed frame to match room decor

4 pc. electro-static mesh bed litter

Electromagnetic powder coated bed frame

4 accessory slots for IV pole, lifting pole or orthopedic frame

Bed exit alarm

CPR function

Night light

Back-up battery (up to 8 hours of battery life)

Central brake

Dual castors

Handheld LINAK nurse control unit

LINAK patient hand controller

Lifting pole

TECHNICAL DETAILS

Size:

(L) 86.6 x (W) 43.3 x (H) 9.8-25.6 (unit of measure: in.)

(L) 2200 x (W) 1100 x (H) 250-650 (unit of measure: mm)

Bed surface area:

(L) 79.9 x (W) 35.4 x (H) 9.8-25.6 (unit of measure: in.)

(L) 2030 x (W) 900 x (H) 250-650 (unit of measure: mm)

Castor size: 5 in. (125 mm) dual castors

Product weight: 160 kg. 300 kg. (661 lbs.)

Safe working load: 200 kg. (440 lbs.)

Back angle: 0-75° (±5°)

Leg angle: 0-35° (±5°)

Trendelenburg & Reverse Trendelenburg angle: 0-13°

ACCESSORIES



